



EN415/EO425 ENG Preamplifier

ENG for balance assessment

Testing flexibility

Interacoustics' Electronystagmography (ENG) option, the EN415/EO425, offers you maximum flexibility when it comes to vestibular testing.

While most patients can be easily tested with the VN415 or VO425 video goggles, you will encounter cases from time to time that require the use of ENG, such as in cases of ptosis (drooping eyelids), or when testing very young patients and this is where EN415/EO425 is useful.

Using electrodes that connect to the EN415/EO425 preamplifier through the computer's USB port, you are able to conduct the relevant tests using the same software you use with your VNG system.

Highlights:

- **Compatible with your VN415/VO425 system**
- **2 or 3 channel (mono or bino) setup**
- **Automatic centering**
- **Easy toggle between ENG and VNG**



Test battery

If you have a VN415, the same test battery will be available with EN415, and if you have a VO425, the same test battery will be available with EO425:

EN415/EO425 and VN415/VO425 tests:

- Bithermal/Monothermal Caloric
- Spontaneous Nystagmus
- Positional
- Dix-Hallpike

EO425 / VO425 include the following additional tests:

- Gaze
- Smooth Pursuit (Tracking)
- Saccade
- Optokinetic



Interacoustics



Interacoustics a/s

Interacoustics is a world leading diagnostic solutions provider in the field of hearing and balance assessment. We help the professional audiological world reach new milestones through continuous developments and a constant focus on integration and direct customer value.

With you at all times

Whenever you purchase audiological equipment with the Interacoustics brand on it, you are not just guaranteed a quality product but also highly trained support personnel with specialized knowledge. We operate in over 100 countries worldwide through a carefully selected network of distributors and service centres to ensure that you receive total support and service.

Product specifications

All technical and hardware specifications concerning the ENG PreAmplifier can be downloaded from our website.

Read more at
interacoustics.com

Go online to
explore our
full product
range

Interacoustics a/s

Drejervaenget 8
DK-5610 Assens

T +45 6371 3555
F +45 6371 3522

info@interacoustics.com
interacoustics.com



3-Channel left binocular setup

